# IT28

Rev. 1 / 12.2020

## Emulator Indicator

IT28, IT28-AD

### Results Reference Guide

**Unprocessed**

**Failures**

One or more parameters of the sterilization cycle were not achieved. The materials cannot be considered sterile and must be processed again.

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
Type 6 Indicator according to ISO 11140-1:2014

### Product Description

Integron® Emulating Indicators have been designed to react to steam sterilization processes inside loads. The blue ink was designed to turn to black when the process reaches the stated values (SV) of the critical process variables indicated in the emulator.

### Characteristics

IT28 consists of a one-side laminated paper strip of 60 mm x 25 mm. IT28-AD consists of a one-side laminated paper strip of 60 mm x 25 mm on a self-adhesive support for record.

Indicators 100 % free of toxic heavy metals.

### Indications for Use

Use the emulator indicator for monitoring Steam sterilization processes. The endpoint is reached in steam sterilization processes of 3.5 minutes at 134 °C or 15 minutes at 121 °C

### Precaution

Do not use this emulator indicator for monitoring Ethylene Oxide, Dry Heat or other sterilization processes different to Steam.

### Instructions of Use

1. Place an Integron® Indicator Emulator within each package, pouch or tray to be sterilized by steam.
2. In rigid containers place an Integron® Emulator Indicator in every corner, or at least two diagonally or opposite corners, to ensure that the steam has reached the most inaccessible areas.
3. After completing the process, analyze the results. If the emulation condition has been reached, the indicator ink must turn to black and it is concluded the sterilization process was successful.

### Storage

Indicators should be stored in dry places, away from light, at temperatures between 10-30 °C, 30-80 % relative humidity. Do not wet. Do not store close to sterilizing agents.

### Shelf Life

Indicators have an expiration date of 5 years from the date of manufacture when stored in the recommended conditions. Both expiration date as well as the batch of the product are specified on the package label and on each unit of product. Do not use indicators after their expiration date.

End Point Stability Reaction: chemical indicator endpoint shall remain unchanged for a period 6 months when stored at previously indicated conditions.

### Disposal

Discard along with cardboard waste according to your country's healthcare regulations.

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Terragene S.A