

Washbon Second Molar Bands

Brand:

Product Code: OWSMB

Availability: Out Of Stock



Price: \$4.20

Ex Tax: \$3.82

10 or more \$3.90

Description

For the ultimate in custom-fitting anatomy, shape and retention.

- Occlusal festoon at distal portion of upper 2nd molar band helps eliminate metal overhang and occlusal interface
- Right and left anatomy is proportioned in all dimensions for positive fit
- Material thickness is tapered from occlusal to gingival, permitting band to pass the contact points when driven

Band No	Upper Left	Upper Right		Band No	Lower Right	Lower Left
97	372-0097	371-0097		1	376-0001	375-0001
98	372-0098	371-0098		2	376-0002	375-0002
99	372-0099	371-0099		3	376-0003	375-0003
00	372-0000	371-0000		4	376-0004	375-0004
1	372-0001	371-0001		5	376-0005	375-0005
2	372-0002	371-0002		6	376-0006	375-0006
3	372-0003	371-0003		7	376-0007	375-0007
4	372-0004	371-0004		8	376-0008	375-0008
5	372-0005	371-0005		9	376-0009	375-0009
6	372-0006	371-0006		10	376-0010	375-0010
7	372-0007	371-0007		11	376-0011	375-0011
8	372-0008	371-0008		12	376-0012	375-0012
9	372-0009	371-0009		13	376-0013	375-0013

10	372-0010	371-0010		14	376-0014	375-0014
11	372-0011	371-0011		15	376-0015	375-0015
12	372-0012	371-0012		16	376-0016	375-0016
13	372-0013	371-0013		17	376-0017	375-0017
14	372-0014	371-0014		18	376-0018	375-0018
15	372-0015	371-0015		19	376-0019	375-0019
16	372-0016	371-0016		20	376-0020	375-0020
17	372-0017	371-0017		21	376-0021	375-0021
18	372-0018	371-0018		22	376-0022	375-0022
19	372-0019	371-0019		23	376-0023	375-0023
20	372-0020	371-0020		24	376-0024	375-0024
21	372-0021	371-0021		25	376-0025	375-0025
22	372-0022	371-0022		26	376-0026	375-0026
23	372-0023	371-0023		27	376-0027	375-0027
24	372-0024	371-0024		28	376-0028	375-0028
25	372-0025	371-0025		29	376-0029	375-0029
26	372-0026	371-0026		30	376-0030	375-0030

Product Gallery

