

NEW YEAR SPECIAL

CATALOGUE

1ST FEBRUARY TO 30TH APRIL, 2026




20/15 Valediction Road Kings Park, NSW 2148

Tel: (02) 9831 6065

hello@dentalzone.com.au | www.dentalzone.com.au

HIGH END COLLAGENS AT AN AFFORDABLE PRICE



STERILE

Indications for use

- Local haemostatis, defect filling into clean wound
- Treatment of extraction socket without bony defect - 4 walled

The 4-in-1 Solution by RESORBA®

- 1. Safe and rapid haemostasis**
 - No risk of secondary bleeding
 - Stabilisation of the blood clot
- 2. Reduction of alveolar ridge atrophy**
 - Stabilisation of the socket and maintenance of the sensitive vestibular lamella
 - Volume stability of the alveolar process
 - Improvement of the aesthetic and functional long-term results
- 3. Reliable bone regeneration**
 - Defined regeneration matrix
 - Angioconductive and osteoconductive properties
 - Already after a relatively short time, a stable bony implant site is formed
 - More easily predictable long-term success
- 4. High user and patient satisfaction**
 - Simple and safe technique
 - Subsequent treatment steps are simplified
 - Higher patient comfort
 - Excellent value for money



Illustrative example of socket healing patterns after extraction.

Above: Socket treated with RESORBA® dental cone.

Below: Untreated alveolar socket.



\$280.00

BOX 10

DK1010
AMS Resorba Parasorb
Collagen Cones
12mm Height



HIGH END COLLAGENS AT AFFORDABLE PRICING



STERILE

Product characteristics of RESORBA® Fleece

- Able to absorb large amounts of fluids on account of their sponge-like structure
- Rapid and safe haemostasis
- Fast integration into the surrounding tissue
- Easy to cut and to adapt

RESORBA® Fleece

- 1 cm² contains 2.8 mg native equine collagen fibrils

RESORBA® Fleece HD

- Highly concentrated collagen for larger bone defects 1 cm² Collagen sponge contains 5.6 mg native equine collagen fibrils



\$333.00

BOX 5

DK9001
AMS Resorba
Parasorb Fleece
70x30 mm

\$387.00

BOX 12

DK1836
AMS Resorba
Parasorb Fleece
18x36 mm

\$440.00


BOX 5

DK8001
AMS Resorba
Parasorb Fleece HD
70x30 mm

\$495.00

BOX 12


DK8002
AMS Resorba
Parasorb Fleece HD
18x36 mm




Absorb-able wound cover Illustrative examples.

Above: RESORBA® Fleece used as a defect filler.

Below: RESORBA® Fleece treating bone void after cyst removal.



COLLAGEN FIXING SUTURES



Needle	Gauge No. USP	Length (Metre)	Colour	Box Content	Code
▼ DS 18	4-0	0.45	Violet	Box - 24	52160
	5-0	0.45	Violet	Box - 24	52507

▼ DS = 3/8 circle, reverse cutting
*Sutures to Scale

AMS Resorba Glycolon® Dental Sutures

Short term absorbable monofilament

- Absorbable mono-filament suture
- Short-term absorption

\$279.00

BOX 24

52160
Glycolon 4/0 18 mm
3/8 Circ R/Cut Violet
45mm DS18

\$317.00

BOX 24

52507
Glycolon 5/0 18 mm
3/8 Circ R/Cut Violet
45mm DS18

COLLAGEN TOOLS



\$15.30

BOX 10

69111
Scalpels #11 Sterile for adjusting
Parasorb Cones



\$17.40

EACH

27120
College Tweezer



\$65.00

EACH

27599
Metzenbaum Fine Scissors
150mm Straight



\$33.00

EACH

27150N
Endo Probe S/E UNC Colour
coded in mm 1-15

NON-LATEX RUBBER DAMS - NEW BULK PACK



NEW BULK PACK

MGuard Non-Latex Rubber Dams - Light Blue

Engineered with enhanced elasticity and greater tear resistance than ever before, the new MGuard Rubber Dams provide superior strength, flexibility, and confidence to every procedure.



\$36.50

68100
MGuard Non Latex Dental Dam
15x15cm Light Blue

EACH

BULK PACK

\$159.50

68100 - 100
MGuard Non Latex Dental Dam
15x15cm Light Blue


100 PACK

duralast

The Single Use Instrument Brand of Choice

MGuard Rubber Dam Frames

DURALAST Rubber Dam Frames are available in both child and adult sizes. The lightweight construction ensures better visibility and easier access during procedures, while the smooth, rounded edges help create a more pleasant experience for patients.



\$17.60

27585
Duralast Rubber Dam Frame Adult

EACH

\$17.60

27586
Duralast Rubber Dam Frame Child

EACH

duralast

The Single Use Instrument Brand of Choice

MGuard Rubber Dam Forceps and Punches

Purchase 1 of the selected Forceps or Punches and receive 4 FREE CLAMPS!

4 free clamp selection (W14A, W2A, W8, 12A, 13A, 26, 2A)

4xFREE CLAMPS WITH PURCHASE*



\$165.00

EACH

27580
Duralast Punch Ainsworth

\$137.50

EACH

27587
Duralast Rubber Dam Forcep Brewer

\$157.00

EACH

27581
Duralast Punch Ivory


\$137.50

EACH

27588
Duralast Rubber Dam Forcep Ivory

AUSTRALIA'S MOST COMPLETE RANGE OF DENTAL NEEDLES

HOGEN **spitze**



\$14.00	\$12.00	\$11.50
BOX 100	BOX 100	BOX 100
EACH	6-11 BOXES	12+ BOXES

- Z1015**
Needles Anaes
Hogenspitze 27G S
21mm Imp Thread
- Z1023**
Needles Anaes
Hogenspitze 27G S
25mm Imp thread
- Z1021**
Needles Anaes
Hogenspitze 27G L
35mm Imp Thread
- Z1019**
Needles Anaes
Hogenspitze 27G L
40mm Imp Thread
- Z1020**
Needles Anaes
Hogenspitze 30G L
25mm Imp Thread
- Z1016**
Needles Anaes
Hogenspitze 30G S
21mm Imp Thread
- Z1022**
Needles Anaes
Hogenspitze 30G XS
12mm Imp Thread

ROUND STEEL BURS

duralast
The Single Use Instrument Brand of Choice

Conveniently sealed in single-cell clean packaging.



- Made from the highest quality alloy steel
- Perfectly balanced to reduce vibration
- Sharp cutting blades
- Economically priced
- Conforms to ISO 3823-1:1997 steel burs standard
- 25 individually sealed burs per box
- ISO CE certified
- TGA listed
- Non-Sterile

\$25.80
BOX 25

ROUND RA Bur diagrams to scale								
Shank Ø	008	010	012	014	016	018	021	023
Round #	#1	#2	#3	#4	#5	#6	#7	#8
SKU	61008	61010	61012	61014	61016	61018	61021	61023

STERILISATION

SELF-SEAL STERI-POUCHES

IN-BUILT INTERNAL & EXTERNAL CLASS 4 INDICATORS

Gives superior sterility information on every pack

No printing within the sealing area No resistance to steam penetration or ink contamination

Medical Grade Paper
As per ISO 11607-1

Wide Triple Seals

More secure construction, less chance of instrument breaching

Blue Tinted Laminate

Identifies any tears or punctures in the film

Score-Lined Ensures accurate, fast seal



\$5.95

BOX 200

1057100M

1. Steri Pouches
Self Seal 57x100

\$8.95

BOX 200

1070230M

2. Steri Pouches
Self Seal 70x230

\$8.75

BOX 200

1090135M

3. Steri Pouches
Self Seal 90x135

\$10.95

BOX 200

1090230M

4. Steri Pouches
Self Seal 90x230

\$17.75

BOX 200

1135260M

5. Steri Pouches
Self Seal 135x260

\$31.50

BOX 200

1190330M

6. Steri Pouches
Self Seal 190x330

\$41.95

BOX 200

1250380M

7. Steri Pouches
Self Seal 250x380

STERILISATION TESTING



\$44.00

PKT-250

54080

SafeCheck Class 6
134°C 3.5 min

- Safe Check Class 6 Sterilisation Control Strips. Conforms to standards ISO 11140-1
- Designed for use in pre-vacuum steam sterilisers with cycle 134°C for 18 minutes
- IFU available.
- TGA Listed for compliance with Australian safety standards.

MEETS AUSTRALIAN GUIDELINE FOR INFECTION PROTECTION & CONTROL



\$59.50

BOX-10

54090

SafeCheck Steam Monitor
Bowie & Dick Test Pack Disposable

- Designed for daily testing of vacuum-assisted steam sterilizers
- Detects air leaks, inadequate air removal, and steam penetration issues
- Easy-to-read colour change indicator for clear pass/fail results
- Single-use, disposable pack for convenience and compliance with standards

MEETS AUSTRALIAN GUIDELINE FOR INFECTION PROTECTION & CONTROL

\$205

KIT

54085

Steam Monitor
SafeCheck Helix Test Kit 1 x PCD + 250 Strips

- Meets EN ISO 11140-6:2022, the standard for testing hollow and porous load sterilisation.
- Designed for pre-vacuum Class B sterilisers to verify steam penetration.
- Certified for 250 cycles, offering reliability and cost efficiency.
- Includes Instructions for Use (IFU) for correct operation and results.
- TGA Listed for compliance with Australian safety standards.

HIGH PERFORMANCE

DENTAGENIE WATER FLOSSER

\$65.00

EACH

GO2WATER
Dentagenie
Water Flosser

A PORTABLE POWERFUL CLEAN

WATERPROOF

At the sink and in the shower too!

MICRO-BUBBLE BLAST

Flushes plaque and bacteria from between teeth and along the gum-line.

CORDLESS

Recharges with USB cable and lasts for up to 1 month.

EASY TO CONTROL

Handy stop/pause button to control water flow easily (once started).

CONSTANT WATER PRESSURE

Even when the battery is low.

SLEEK DESIGN

Modeled around the shape of your hand and looks fabulous on your vanity.

SMART TIMER

Automatically powers off when the water tank is empty or 2 minutes is up.



INCLUDES:

- 1 WATER FLOSSER
- 1 NOZZLE
- 1 USB CHARGING CORD

3 FLOSSING MODES




SOFT MODE

Gently cleans around
braces... Easy!



CLEAN MODE

Blasts between teeth...
Super clean!



PULSE MODE

Cleans and massages
gums... Love it!



PERIOGOLD

Helps Improve Oral Health

With Microencapsulated Curcumin, Aloe Vera & Elderberry

TIME TO REPLACE CHLORHEXIDINE WITH PERIOGOLD

A safer, gentler, and more effective approach to modern oral care.

A Natural, Herbal-Based Solution for Healthier Gums and Happier Patients



Benefits of using PerioGold over Chlorhexidine

- **No staining or taste disturbance** – avoids the common patient complaints associated with Chlorhexidine.
- **Gentle on oral tissues** – reduces irritation and mucosal damage often linked to prolonged Chlorhexidine use.
- **Herbal formulation** – powered by scientifically supported botanicals traditionally used for antimicrobial and anti-inflammatory benefits.
- **Effective antimicrobial action** – targets harmful oral pathogens without disrupting overall oral flora balance.
- **Better patient compliance** – pleasant flavour and non-harsh feel mean patients are more likely to use it correctly and consistently.
- **Safe for longer-term use** – suitable for extended periodontal maintenance regimes unlike Chlorhexidine's limited-use recommendation.
- **Supports holistic care** – aligns with the growing patient demand for natural or plant-based oral care options.

10075 PerioGold Gel 40mL (46g)

\$17.70

EACH

\$15.40

6-11 BOXES

\$13.70

12+ BOXES

Floss



Brush



Restore

