

1ST MAY TO 30TH JUNE, 2026

EOFY SALE

CATALOGUE



20/15 Valediction Road Kings Park, NSW 2148

Tel: (02) 9831 6065

hello@dentalzone.com.au | www.dentalzone.com.au



STERILE

FROM SMART SYNTHETICS TO NATURAL BONE

NEW
PRODUCT

RESORBA® Synthetic Bone

YOUR VIABLE OPTION



NEW

Biphasic Formulation
B - Tricalcium Phosphate (B-TCP)
+ Hydroxyapatite (H)



\$198

1 SYRINGE

SB0401GS50
RESORBA® Synthetic Bone –
Granules in a syringe 60%
HA/40% β -TCP mix
Particle Size: 0.5 mm – 1.0 mm
Box Contents: 1 syringe, 0.5 cc

Key Features

- Maximum permeability
- 70% global porosity – similar to cancellous bone
- Easy to use
- Different shapes & ratio available
- Osteoconductive – provides a scaffold for new bone growth¹
- Biphasic composition of HA and TCP provides optimum performances for a gradual replacement of the substitute by human bone

- ✓ 30+ years of clinical experience
- ✓ More than 1 million units manufactured

Therapeutic Suggestions

- Use as an autograft extender
- Implantology:
 - Guided Bone Regeneration procedures (GBR)
 - Sinus-lift augmentation
 - Vertical ridge augmentation
 - Horizontal alveolar ridge augmentation
 - Alveolar regeneration – Extraction socket
- Periodontology:
 - Intra-osseus pockets
 - Furcations


1. Rodriguez C, Jean A, Mitja S, Daculsi G. Five Years Clinical Follow-Up. Bone Regeneration with CaP Bioceramics. Key Engineering Materials. 2008;361-363:1339-42.



STERILE

NEW
PRODUCT

FLEXIBLE COLLAGEN MEMBRANE WITH CONTROLLED RESORPTION TIME



- Adaptable** – dry or wet
- Strong** – easy to suture
- Cross-linked** – for desired barrier function
- Cuttable** – size or shape can be adjusted with scissors

Easy to use
on both sides

Product Description

RESORBA® EasyFlex is a resorbable membrane made from porcine dermis-sourced collagen, aimed for guided bone and tissue regeneration (GBR & GTR).

Due to its substantially proven biocompatibility and low immunogenicity, porcine has been

used in a variety of maxillofacial procedures for many years.

RESORBA® EasyFlex is a cross-linked collagen membrane that combines a predictable resorption time and flexibility. Providing the desirable barrier durability for 12 weeks.


Key Indications

RESORBA® EasyFlex is a resorbable membrane. It is intended for use in periodontal/dental surgery procedures, for placement in areas of periodontal defect, dental implant, bone defect or ridge reconstruction; to support wound healing post surgery.

RESORBA® EasyFlex is intended to create a physical barrier at the surgical site to prevent the

colonisation of the underlying bone site by soft tissue cells (connective and epithelial) for 12 weeks while allowing the diffusion of biological fluids. As the membrane is naturally resorbable, there is no need for a second surgical intervention to remove it. The total resorption time is 24 weeks.

RESORBA® EasyFlex is recommended for use with a bone graft material to promote new bone healing.



\$269

1 BOX

SB0701EZC1525

RESORBA® EasyFlex
PORCINE CO

Size: 15 x 25 mm (3,75 cm²)

Contents: 1 membrane



NEW PRODUCT

HIGH END COLLAGENS AT AN AFFORDABLE PRICE



STERILE

Indications for use

- Local haemostatis, defect filling into clean wound
- Treatment of extraction socket without bony defect - 4 walled

The 4-in-1 Solution by RESORBA®

1. Safe and rapid haemostasis
 - No risk of secondary bleeding
 - Stabilisation of the blood clot
2. Reduction of alveolar ridge atrophy
 - Stabilisation of the socket and maintenance of the sensitive vestibular lamella
 - Volume stability of the alveolar process
 - Improvement of the aesthetic and functional long-term results
3. Reliable bone regeneration
 - Defined regeneration matrix
 - Angioconductive and osteoconductive properties
 - Already after a relatively short time, a stable bony implant site is formed
 - More easily predictable long-term success
4. High user and patient satisfaction
 - Simple and safe technique
 - Subsequent treatment steps are simplified
 - Higher patient comfort
 - Excellent value for money



Illustrative example of socket healing patterns after extraction.

Above: Socket treated with RESORBA® dental cone.

Below: Untreated alveolar socket.



\$259
 BOX 10

DK1010
 AMS Resorba® Parasorb
 Collagen Cones
 Size: 12mm Height
 Contents: 10 Cones



ECO-FRIENDLY BARRIER PROTECTION FOR EVERY SURGERY




\$14.00
BOX 500

85071E
Barrier Cover **BIO** Syringe Sleeve with Opening



\$15.90
BOX 500

81003E
Barrier Cover **BIO** H/Speed Handpiece - 25 x 203mm



\$14.00
BOX 500

81004E
Barrier Cover **BIO** S/Speed Handpiece 38mm x 203mm



\$28.00
BOX 400

85093E
Barrier Cover **BIO** LED Curing Light 330 x 53mm



\$28.00
BOX 250

85084E
Barrier Cover **BIO** Headrest 250 x 280mm



\$28.00
BOX 250

85085E
Barrier Cover **BIO** Headrest Covers 254 x 355mm



\$8.90
BOX 200

85096E
Barrier Cover **BIO** Light Guide Tip 7-8mm



\$8.90
BOX 200

85097E
Barrier Cover **BIO** Light Guide Tip 9-11mm



\$15.90
BOX 500

85099E
Barrier Cover **BIO** Scaler 40 x 203mm




\$15.90
BOX 500

85098E
Barrier Cover **BIO** Scaler 43 x 177mm



\$38.90
BOX 500

85081E
Barrier Cover **BIO** Tray 200 x 270mm F Tray



\$41.50
BOX 300

85080E
Barrier Cover **BIO** Tray Sleeves 266 x 355mm




\$11.50
BOX 500

85091E
Barrier Cover **BIO** X-ray Sensor 41 x 203mm



\$13.90
BOX 500

85092E
Barrier Cover **BIO** X-ray Sensor 60 x 200mm



\$12.90
BOX 500

85090E
Barrier Cover **BIO** T Bar Light Handle



\$38.90
BOX 250

85075E
Barrier Cover **BIO** X-Ray Head 38cm x 66cm



\$28.00
BOX 100

80097
Barrier Cover Intra Oral Camera Sopro 717

Barrier Covers

- Defense from contamination provided by the envelope allows PSP sensors to be reused.
- Black opaque envelope back prevents ambient light exposure resulting in distorted images.

\$29.70
BOX 200

53060	PSP Size 0 Short Edge Loading
53999	PSP Size 1 Long Edge Loading
54001	PSP Size 1 Short Edge Loading
54000	PSP Size 2 Long Edge Loading
54008	PSP Size 2 Short Edge Loading
53063	PSP Size 3 Short Edge Loading
53064	PSP Size 4 Short Edge Loading

Barrier Film Non-Stick Edge



\$20.00
ROLL

81000 Barrier Film - Blue
81001 Barrier Film - Clear



\$89.90

CTN 2000

200204
 Cotton Rolls Braided 10mm #2



\$17.00

BOX 1000

39014
 Cotton Rolls Celleron
 #1 8mm



\$15.00

BOX 1000

39015
 Cotton Rolls Celleron
 #2 10mm



**BUY
 5 GET
 1 FREE**

\$39.90

BOX 50

49001
 Syringe Endo Irrigating/
 Utility Curved Tip



**BUY
 5 GET
 1 FREE**

\$32.00

BAG 250

49210 Clear Syringe Air Water
 49211 White Disposable Tip Accu Tip
 49212 Blue



**BUY
 5 GET
 1 FREE**

\$32.00

BAG 100

49000
 Syringe Air Water Disposable
 Tip Denti-Tips White



**BUY
 5 GET
 1 FREE**

\$32.00

BAG 250

49200
 Syringe Air Water Disposable
 Tip Accu Tip - Colour

SELF-SEAL STERI-POUCHES



IN-BUILT INTERNAL & EXTERNAL CLASS 4 INDICATORS
Gives superior sterility information on every pack


No printing within the sealing area No resistance to steam penetration or ink contamination

Medical Grade Paper
As per ISO 11607-1

Wide Triple Seals
More secure construction, less chance of instrument breaching

Blue Tinted Laminate
Identifies any tears or punctures in the film

Score-Lined Ensures accurate, fast seal



\$5.90 BOX 200	\$8.95 BOX 200	\$7.90 BOX 200	\$9.99 BOX 200	\$15.50 BOX 200	\$27.50 BOX 200	\$38.90 BOX 200
1057100M 1. Steri Pouches Self Seal 57x100	1070230M 2. Steri Pouches Self Seal 70x230	1090135M 3. Steri Pouches Self Seal 90x135	1090230M 4. Steri Pouches Self Seal 90x230	1135260M 5. Steri Pouches Self Seal 135x260	1190330M 6. Steri Pouches Self Seal 190x330	1250380M 7. Steri Pouches Self Seal 250x380

STERILISATION TESTING



\$46.00 1+ PKTS
PKT 250

\$42.00 3+ PKTS
PKT 250

MEETS AUSTRALIAN GUIDELINE FOR INFECTION PROTECTION & CONTROL

54080
SafeCheck Class 6
134°C 3.5 min

- Safe Check Class 6 Sterilisation Control Strips. Conforms to standards ISO 11140-1
- Designed for use in pre-vacuum steam sterilisers with cycle 134°C for 3.5 minutes
- IFU available.
- TGA Listed for compliance with Australian safety standards.



\$48.00
BOX 10

MEETS AUSTRALIAN GUIDELINE FOR INFECTION PROTECTION & CONTROL

54090
SafeCheck Steam Monitor
Bowie & Dick Test Pack Disposable

- Designed for daily testing of vacuum-assisted steam sterilizers
- Detects air leaks, inadequate air removal, and steam penetration issues
- Easy-to-read colour change indicator for clear pass/fail results
- Single-use, disposable pack for convenience and compliance with standards




\$179
KIT

54085
Steam Monitor
SafeCheck Helix Test Kit 1 x PCD + 250 Strips

- Meets EN ISO 11140-6:2022, the standard for testing hollow and porous load sterilisation.
- Designed for pre-vacuum Class B sterilisers to verify steam penetration.
- Certified for 250 cycles, offering reliability and cost efficiency.
- Includes Instructions for Use (IFU) for correct operation and results.
- TGA Listed for compliance with Australian safety standards.

AUSTRALIA'S MOST COMPLETE RANGE OF DENTAL NEEDLES

HOGENspitze



\$11.90
 BOX 100
 1 - 9 BOXES

\$11.50
 BOX 100
 10+ MIX BOXES

- Z1015**
Needles Anaes
Hogenspitze 27G S
21mm Imp Thread
- Z1023**
Needles Anaes
Hogenspitze 27G S
25mm Imp thread
- Z1021**
Needles Anaes
Hogenspitze 27G L
35mm Imp Thread
- Z1019**
Needles Anaes
Hogenspitze 27G L
40mm Imp Thread
- Z1022**
Needles Anaes
Hogenspitze 30G XS
12mm Imp Thread
- Z1016**
Needles Anaes
Hogenspitze 30G S
21mm Imp Thread
- Z1020**
Needles Anaes
Hogenspitze 30G L
25mm Imp Thread

ROUND STEEL BURS

Conveniently sealed in single-cell clean packaging.

duralast
 The Single Use Instrument Brand of Choice



- Made from the highest quality alloy steel
- Perfectly balanced to reduce vibration
- Sharp cutting blades
- Economically priced
- Conforms to ISO 3823-1:1997 steel burs standard
- 25 individually sealed burs per box
- ISO CE certified
- TGA listed
- Non-Sterile



\$28.70
 BOX 25


ROUND RA Bur diagrams to scale								
Shank Ø	008	010	012	014	016	018	021	023
Round #	#1	#2	#3	#4	#5	#6	#7	#8
SKU	61008	61010	61012	61014	61016	61018	61021	61023

4+1 ULTRA ABSORBENT BIBS



Deliver comfort and protection with our top-quality MGuard bibs — designed for performance and convenience.

- **5-PLY PROTECTION:**
4-ply soft paper + 1-ply polyethylene
- **ULTRA-ABSORBENT**, soft texture for maximum comfort
- **EMBOSSED** for superior fluid retention
- **REINFORCED EDGES** for added strength
- Packed in a sealed, easy-access drawer pouches (100 sheets)
- Available in small and large sizes
- Perfect for use as patient bibs or tray liners



CONVENIENT DRAWER POUCHES CONTAINING 100 BIBS PER POUCH



\$49.00
BOX 800

40107 Small Bib MGuard, 4-Ply Paper + 1 Ply Polyethylene, White, 28.5 cm x 21.5 cm



\$49.00
BOX 400

40106 Large Bib, 4-Ply Paper + 1 Ply Polyethylene, White, 28.5 cm x 43.0 cm

FACIAL AND BODY PROTECTION



Level 2 Barrier - 120mmHg

- Exceeds Australian Standard AS 4381.2015
- Level 2 Barrier Protection: For low to moderate levels of aerosols, spray and/or fluids
- BFE Bacteria Filtration Efficiency 3 micron > 95%
- PFE Particle Filtration Efficiency 0.1 micron > 98%
- Listed on the TGA

Anti-Fog Level 2 Barrier - 120mmHg

- Exceeds Australian Standard AS 4381.2015
- Level 2 Barrier Protection: For low to moderate levels of aerosols, spray and/or fluids
- BFE Bacteria Filtration Efficiency 3 microns > 95%
- Anti-Fog Filter
- PFE Particle Filtration Efficiency 0.1 micron > 98%
- Listed on TGA

Anti-Fog Level 3 Barrier - 160mmHg

- Exceeds Australian Standard AS 4381.2015
- Level 3 Barrier Protection: For high risk of fluid, spray or aerosol exposure
- BFE Bacteria Filtration Efficiency 3 microns > 95%
- Anti-Fog Filter
- PFE Particle Filtration Efficiency 0.1 micron > 98%
- Listed on TGA

PRODUCT CODE	FACE MASK BARRIER LEVEL	FACE MASK TYPE	ANTI-FOG	MASKS PER BOX	PRICE 1-19 UNITS	PRICE 20+ UNITS
70200	2	Tie on		50	\$6.00	\$6.00
70204	2	Earloop		50	\$4.50	\$4.00
70206	2	Tie on	✓	50	\$6.90	\$6.50
70208	2	Earloop	✓	50	\$5.90	\$5.50
70214	3	Earloop	✓	50	\$12.50	\$11.90



\$22.00
10 PACK

70110
Gown MGuard SMS Blue AAMI Level 1 Protection

- Infection Control gown with neck and waist ties and cotton stretch cuff
- Listed on the TGA (334546)



\$40.00
PKT 10

70218
Face Shield w Anti-Fog Coating

- Full Length Face Shield
- Anti-Fog, Anti-Static and Anti-Glare
- Listed on the TGA (331423)